

## SPEC SHEET

# Inspection Table

## JYT-A MODEL

### INTRODUCTION

Allows for effective inspection and selection, for the removal of unqualified capsules to ensure they are not mixed with the finished product.

**KEYWORDS** [Inspection table](#) | [inspection and selection](#) | [unqualified softgel capsule](#) | [softgel](#) | [accessories](#)



### ATTRIBUTES

- Waste pill bucket for singling out unqualified capsules
- Lighting lamp for clear inspection

### TECHNICAL PARAMETERS

Model	JYT-A
Power Supply	220V 60w
Dimensions	1200×550×1200mm
Weight	65kg

While the information has been carefully compiled to the best of our knowledge, nothing is intended as representation or warranty on our part and no statement herein shall be construed as recommendation to infringe existing patents. In the effort to improve our products, we reserve the right to make changes judged to be necessary.

### CONNECT WITH US

+86 577 6688 6278 | [info@capsulatech.com](mailto:info@capsulatech.com) | [www.captksolutions.com](http://www.captksolutions.com)